
Resonant Frequency Doubler for CW single-frequency lasers, model "FD-SF-07"

Tekhnoscan presents **resonant frequency doubler**, **model FD-SF-07**, with **Smart Auto-Relock function** for CW single-frequency lasers (solid-state, fiber, dye, etc.) that opens a new possibilities for more efficient laser wavelength conversion in the visible and near IR ranges into the blue and UV domains. **Optimised resonator** of FD-SF-07 in combination with **high-quality mirrors** ensures relatively high level of output second-harmonic power. Pumped with **1 W** fundamental radiation power the doubler outputs: **more than 250 mW** within the 700-950 nm range, **> 200 mW** within the 550-700 nm range, and **> 150 mW** within the 550-500 nm range.



The Smart Auto-Relock function allows FD-SF-07 doubler to smoothly track considerable changes in the frequency of the input radiation, thus the range of smooth second-harmonic frequency scanning may cover **dozens of GHz**, being only limited by the spectral acceptance of the non-linear crystal. The FD-SF-07 doubler is notable for its low acoustic noise and sensitivity to vibrations, as well as for the simplicity of tuning and ease of use. **Super-stable and compact ring cavity** combined with ultra-fast two-stage system that locks the cavity to the frequency of the input radiation by the Hansch-Couillaud method are a guarantee for **high stability of the output power** of the second harmonics even for lasers without a frequency stabilisation.



Advanced Realized
Photonics Ideas



- Improved doubling efficiency up to 40% at the input radiation power 1 W
- Power-enhancement factor up to 130
- Rigid and stable, easily aligned cavity with vibration isolation



- Fast two-stage system of locking the cavity to the frequency of input radiation
- Possibility of efficient operation with lasers without frequency stabilisation
- Solid ultra-stable performance even under conditions of considerable external acoustic perturbations and vibrations
- Thermo-stabilised non-linear crystal



- A separate control photodetector included into the standard set of the electronic unit
- © Convenient express-access to the alignment of the nonlinear crystal through a special easily removable lid

FD-SF-07

Resonant Frequency Doubler

Specifications:

Conversion efficiency for 1 W CW single-frequency input:

700-950 nm: > 25% 550-700 nm: > 20% 500-550 nm: > 15%

Nonlinear crystal types: LBO / BBO

Typical tuning range


for one optics set: 45-75 nm for fundamental input

Features:


- unique function Smart Auto-Relock
- ▶ temperature stabilization of crystal
- ▶ input coupling optics
- output beam-shaping optics


Options:

1. Brewster-angle crystal design for definite wavelength (532 nm, etc.)




Typical input/output curves






Dimensions:














Information and specifications contained herein are deemed to be reliable and accurate as of the publication date. Tekhnoscan reserves the right to change these specifications at any time without notice.



Company Headquarters Sirenevaia Str., 37, k. 141, Novosibirsk, 630058 Russia Tel.: +7-(383)-214-00-09 Fax: +7-(383)-363-42-65 E-mail: service@tekhnoscan.com

www.tekhnoscan.com



call us, or visit

tekhnoscan.com/english/coordinates.htm





